WashU Epigenome Browser Documentation

Release 1.0

Daofeng Li, Silas Hsu, Renee L. Sears and Ting Wang

Nov 26, 2018

Contents:

1	Use t	he Browser
	1.1	The Top Menu
	1.2	Genomic Region Locator
	1.3	Operations in the Tracks View Container
	1.4	Settings
	1.5	Apps 10
	1.6	Track management
	1.7	Track Customization
	1.8	Circlet view for chromatin interaction tracks
2	Trac	KS 23
	2.1	Configure your webserver to enable CORS 23
	2.2	Prepare track files
	2.3	Annotation Tracks
	2.4	Numerical Tracks
	2.5	Read Alignment BAM Tracks 27
	2.6	Methylation tracks
	2.7	Categorical tracks
	2.8	Long range chromatin interaction
3	Data	hub 31
	3.1	Example data hub
	3.2	Example bigWig track
	3.3	Example methylC track
	3.4	Example categorical track
	3.5	Example longrange track
	3.6	Example bigInteract track
	3.7	Example repeatmasker track
	3.8	Example geneAnnotation track 34
	3.9	Example bigbed track
	3.10	Example bed track
	3.11	Example refbed track
	3.12	Example HiC track
	3.13	Example cool track
	3.14	Example genomealign track
	3.15	Example Ruler track
	3.16	Track properties

4	Installation4.1Setup4.2Start the browser4.3Setup your own backend API (optional)	. 41
5	Embedding	43
6	Frontend code architeture	45
	6.1 Quick tour	. 45
	6.2 Suggested order of reading	
	6.3 Making a new track type	
	6.4 Performance tips	. 47
	6.5 React (and other) gotchas	. 47
	6.6 Lessons trying to refactor into WebWorkers	. 47
7	Contact us	49
	7.1 Source code and issue tracker	. 49
	7.2 Get help on social media	. 49
8	Indices and tables	51




Fig. 1: Gateway to the epigenome. (Art by Ting Wang)

CHAPTER 1


Use the Browser

1.1 The Top Menu

The New v47.2 WASHU EPIGENOME BROWSER	Mouse mm10	chr6:52425276-52425961	+%	-%	-1	-5	Tracks 🕶	🖉 Apps 🔻	ØSettings 🕶	

The top navigation menu (above) controls most of the browser functionality. From left to right are the browser logo with version information, species and assembly information, genomic region locator, zoom in/out tools, Tracks menu, Apps menu, Settings menu and Documentaion link.


1.2 Genomic Region Locator



The genomic region locator allows the user to navigate to a region or gene.

1.2.1 Gene Search

You can type a gene name/symbol, like Hox, and when the input content reaches 3 characters the browser will try to find gene symbols starting with what you typed:



After a gene is selected a dropdown menu will pop up with isoforms for the gene. After clicking an isoform the browser will navigate you to its genomic region.

Gene searc	h	Say a Gene Reset Stop	×		
Hoxc10					
refGene	chr15:102966795- 102971897	• > > > > > > > > > > > > > > > > > > >	Mus mRN	musculus homeobox C10 (IA.	Hoxc10),
gencodeM18	chr15:102966796- 102971893	>>> >>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>	Mus mRN	musculus homeobox C10 (IA.	Hoxc10),

Voice input gene symbol

From this set of buttons, you for permission to access your microphone devices, choose *Allow*, and the browser will start to listen to what you are saying. You can start saying letters one by one, like H, O, X, if you click the red **stop** button, what you said "HOX" will populate the gene search box and suggested gene symbols will pop up. As before you can then choose the gene and isoform you want to navigate to.

Note: This feature is dependent on web browser support. A web browser without support for speech recognition won't see this UI.

1.2.2 Region Search

Below the gene search box you can use the region search box to navigate to specific genomic coordinates. Formats such as chr6:52258852-52260880 or chr6 52258852 52260880 are accepted (the browser is not sensitive to the number of spaces or tabs between the *chr*, *start*, and *stop*.

	Region search (current region is chr6:52258852-52260880)									
-	chr1	0	10000	Go	-					

1.3 Operations in the Tracks View Container



Right above the tracks on the left hand side there are a bunch of tools for operating on the tracks. From left to right these include the move, re-order, zoom, undo, redo, and history tools.

1.3.1 Move/Re-order/Zoom tools

The is the default tool selected by the browser. Moving mode allows the user to drag the mouse right or left and the tracks will move along mouse moving to new regions. The allows tracks to be re-ordered. The user can drag one or more tracks up or down to change the order of tracks. The allows the user to zoom in on a specific region within the current view using the mouse.

1.3.2 Undo/Redo/History tools

contains the Undo, Redo, and History tools. For instance if you accidently moved to another region and forgot what you were looking at before you can click the Undo button to go back. Clicking the Redo button allows you to go forward to the step before you clicked Undo. The History button gives you the 9 most recent operations and allows you to jump to any of these operations or clear the history.

Operation history

Close Clear History

Go back:

- 1. Region: chr6:52425276-52425961, # of tracks: 4
- 2. Region: chr6:52258852-52260880, # of tracks: 4
- 3. Region: chr6:52258852-52260880, # of tracks: 4

Go forward:

- 1. Region: chr6:52257584-52259612, # of tracks: 4
- 2. Region: chr6:52257998-52260026, # of tracks: 4

1.3.3 Hotkeys

- 1. Alt + H or Alt + D for the Drag Tool
- 2. Alt + S or Alt + R for the Reorder/Swap Tool
- 3. Alt + M for the Magnify Tool
- 4. Alt + Z and Alt + X to pan one full panel left or right.

1.4 Settings

The Settings menu controls global settings for the browser.

1.4.1 Toggle display of the Genome Navigator

By using the genome navigator (below) users can jump to any genomic region or chromosome(s).



The operations on the genome navigator are:

- Left mouse drag: select
- Right mouse drag: pan
- Mousewheel: zoom

The genome navigator can also be hidden to save space when viewing tracks. Click Settings on the top menu and uncheck the box to switch off this feature:



1.4.2 Toggle highlighting of enter region

When a user jumps to a region or gene using the Genomic Region Locator, that region or gene is highlighted with a light yellow box.

Ο	C										
-		Hoxa10 Hoxa10				<	<u> </u>		←← ╡ < < < ←← ← ← ← ←	< < < < < < < < < <	< < < < < <
							рВ3				
К		52231K		1	52232K		52233K		52234K	5	2235K
							chr6				
_				•						-	
	< <hox< td=""><td></td><td></td><td></td><td>0< < < < <</td><td>< < < < <</td><td>~~~~</td><td></td><td>****</td><td>~~~~~</td><td>····</td></hox<>				0< < < < <	< < < < <	~~~~		****	~~~~~	····
		Hoxa10			xa10< < < < < < < <		$\overline{\cdots}$	\leftrightarrow	~~~ ~~~	<<<<<	
<_H	oxa9	Hoxa10	K < < •	< < < <	< < < < <	< < < < <	~~~ <	l i			

This highlighting can be turned off/on by clicking the botton on the Settings menu:

♦Settings ▼ □Documentation
□ Show genome-wide navigator?
☑ Highlight entered region

1.4.3 Change track label width

The default width of track labels (below) is 120 pixels.


refGene				
Ruler		1	1	
RepeatMasker	1.0 -			
	0.0 -			
gencodeM18		< < 		< < <
		~		~~

The width of the track label can be configured by the submenu under the Settings menu:


ØSettings ▼								
Show genome-wide navigator?								
🗷 Highlig	ht entered region							
Show 3D scene								
120 🌻 🕻	Change track legend width							

1.4.4 Toggle display of VR mode

From the Settings menu the user can choose to toggle the VR display mode of tracks:




والمالية المراجع والمحدية والمراجعة والمراجعة والمراجع والمحدية والم



After choose the Show 3D scene submenu, a new container with VR view of the tracks will appear:




You can click the source icon at the bottom right to toggle the full screen display of VR mode, then you can use your mouse and keys W, A, S and D to control the view of VR mode, like this view below can show you the interaction between two genomic loci and methylation status along this region in a 3D way.



1.5 Apps

1.5.1 Region set view

Users can submit a list of regions or genes to the browser, by choose Apps -> Region set view:



The brings up the region set user interface, here you can enter a list of gene names or coordinates to make a gene set one item per line. Gene names and coordinates can be mixed for input. Coordinate string must be in the form of "chr1:345-678" fields can be joined by space/tab/comma/colon/hyphen.

Select a gene/region set

Add new set

Create a new set

Enter a list of regions

Enter a list of gene names or coordinates to make a gene set one space/tab/comma/colon/hyphen.

CYP4A22	
chr10:96796528-96829254	
CYP2A6	
CYP3A4	
chr1:47223509-47276522	
CYP1A2	
	//
	//
Add Clear	

After Click the *Add* button, will bring you to the region set editting interface, you can either add region one by one, or delete regions from the table, and set the flanking region strategy:

Select a gene/regi	on set			
Add new set				
Create a new set				
1. Rename this set: New set				
2. Add one region or delete regi	on(s) from the table below			
New region name:	New region locus:	Add new region		
Name	Locus	Strand	Coordinates to view	
chr10:96796528-96829254	chr10:96796528-96829254		chr10:96796528-96829254	Delete
chr1:47223509-47276522	chr1:47223509-47276522	-	chr1:47223509-47276522	Delete
Cyp1a2	chr9:57676936-57683655		chr9:57676936-57683655	Delete
Pre	vious Page	1 of 1	10 rows 🜲	Ne
3. Set flanking region				
Upstream bases: 0	Downstream bases: 0	Surrounding: Gene body		

Once you done with edit the set, cick the button Add set. Now you have the option to enter regin set view, click the button Enter region set view:

Select a gene/region set New set (3 regions) Enter region set view DELETE This indicates you are in *region set view* mode and which set you are viewing: Exit region set view Is current view New set (3 regions)

Go back to the browser, you can your browser view is ordered by your region set:

	47.2.1 OWSER Mouse mm10 chr10:96796528-96829254:0-0	Cyp1a2:6719 +½ -⅓ -1	-5 Tracks - 🔨 Ap	ps 🕶 🔍 Settings 👻 🛄 Doc	cumentation		
96800K	qC3 96810K 96820K	47230K	47240K	qC1.1 47250K	47260K	47270K	qB 57680K
	chr10:96796528-96829254			chr1:47223509-47276522			Cyp1a2
Tools: 🖖 🖾	□9+ 0 0 ■						Metadata »
refGene						Сур	1a2 ⊳→⊠>I>I∜ >I>
gencodeM18	≅RP23-40N7.4	→ → → → → → · RP23-40N7.3				Сур	1a2 ►>∎X>I>I≾→I>→ Cyp1a2® -N
gencodeM18Basic	BRP23-40N7.4	>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>				Сур	1a2 >>>>>1>1\$ > > >

DELETE

1.5.2 Session

Choosing Session from the Apps menu will bring you to the session interface shown below:

Session bundle Id	Retrieve session	×						
Save session Session bundle Id: 9d99b640-c4c8-11e8-8121-0d143821f96b								
Name your session: Incredible-amber-zonkey	or use a	Random name						

Save session

Click the **Save session** button to save a session. A session bundle Id will be created which allows the user to retrieve their session at a later date.

Session bundle Id			esession	×			
Save session Session bundle Id: 9	d99b640-c4c8-11e8-8121-0d1	43821f96	ìb				
Name your session:	Whimsical-peach-lynx		or use a	Random name			
1. Spotted-copper-oyster (9/30/2018, 10:58:38 AM) Restored Delete 2. Incredible-amber-zonkey (9/30/2018, 10:58:24 AM) Restore Delete							

Retrieve session

The **session bundle Id** can be used later to retrieve a session by pasting the session bundle id in the session interface and clicking the Retrieve session button.

9d99b640-c4c8-11e8	8-8121-0d143821f96b	Retrieve session	×
Save session Session bundle Id: e	e71b4e0-c4c9-11e8-9baa-03c	ce76b7482b	
Name your session:	Scientific-bronze-clam	or use a	Random name
6			ľ

Choose which session status you want to restore:


- 1. Spotted-copper-oyster (9/30/2018, 10:58:38 AM) Restore Delete
- 2. Incredible-amber-zonkey (9/30/2018, 10:58:24 AM) Restore Delete

Click the green Restore button and your session will be restored:

٦v	The New v47.2 VASHU EPIGENOME BROWSER Mouse mm10 chr6:52424728-52425413 +% -%	-1
	9d99b640-c4c8-11e8-8121-0d143821f96b Retrieve session ×	
	Save session	10
	Session bundle Id: 9d99b640-c4c8-11e8-8121-0d143821f96b	chr
Тоо	Name your session: Scientific-bronze-clam or use a Random name	
refGe		
geno Ruler	2. Incredible-amber-zonkey (9/30/2018, 10.58.24 Alvi) Restore Delete	
	52424800 52424900	
Repe	eatMasker 1.0	

1.5.3 Live browsing

From the Apps menu choose **Go Live**, the browser will navigate you to a new link which you can share with someone else, like your collaborator, your PI, or your friends. Whatever operations are done by you are mirrored on the displays of the people who opened the same link.



1.5.4 Screenshot

Users can create publication quality images using the *Screenshot* tool from the Apps menu. Click the *Screenshot* button and a new window will po pup that re-renders all your tracks as a new SVG file. Once rendered you can click the green download button to save the current browser view as a SVG image file.

Ruler oM 52.1M 52.2M 52.3M 52.4M 52.5M						chr6	
Ruler		OM	52.1M	52.2M	52.3M		52.5M
* Download	Ruler					qB3	
Description	↓ Download						

1.6 Track management

The browser collects data from large corsortia like Roadmap Epigenomics, ENCODE, 4DN, TaRGET, etc. The data are called public data/tracks and are organized into different collections called hubs. Along with these public hubs and tracks users can submit their own custom tracks and data hubs to allow for easy comparison.

1.6.1 Add tracks from public hubs

From the Tracks menu choose **Public Data Hubs**. This will display all of the public data hubbs available for the species and build you are currently working in. For example, using mouse mm10 annotation the *4D Nucleome Network* hub is available. Click the *Add* button to load this hub:

Public data hubs

Collection	Hub name	Tracks	Add
4D Nucleome Network	4DN HiC datasets	23	+

Collection details

The 4D Nucleome Network aims to understand the principles underlying nuclear organization in space and time, the role nuclear organization plays in gene expression and cellular function, and how changes in nuclear organization affect normal development as well as various diseases. The program is developing novel tools to explore the dynamic nuclear architecture and its role in gene expression programs, models to examine the relationship between nuclear organization and function, and reference maps of nucleararchitecture in a variety of cells and tissues as a community resource.

After a hub is added, a facet table containing all tracks will pop up. This allows you to choose any tracks you are interested in. The facet table can also be revisted through the menu when you choose **Track Facet Table**:

Previous	Page 1 of 1	10 rows \$	Next
Row: Sample \$	ك		Column: Assay
			⊞Assay
⊞Sample			0/23

You can expand the row and/or column selection by clicking the + buttons. Row and column displays can also be easily swapped:


Row: Sample 🔶 🗢	\$		Column: Assay			
	BAssay	EDNase Hi-C	Enzyme: DNasel	Bin situ Hi-C	Enzyme: Mbol	Enzyme: Dpnll
⊟Sample						
⊟tissue						
brain			0 /1			
<i>⊟stem cell</i>						
ES-E14					0 /3	
Stable Transfection complex change for Gene:MLL3,MLL4					0/3	
<i>⊟immortalized cell line</i>						
CH12.LX					0 /1	
Patski			0 /5			
⊟stem cell derived cell line						
F123-CASTx129 (Tier 2)					0 /4	
⊟primary cell						
globose cell of olfactory epithelium						0 /1
immediate neuronal precursor cell						0 /1
olfactory receptor cell						0/4

Clicking a cell within the facet table will pop up a new window containing a table with the tacks that match the row and column selections:

Track table

Name	Data hub	Sample	Assay	Format	Add
		Filter	Filter		
brain DNase Hi-C	4DN HiC datasets	tissue > brain	DNase Hi-C > Enzyme: DN	hic	+

Click the *Add* button to add the track(s) you want. You can then view tracks in the browser view window:



×

1.6.2 Adding annotation tracks

Users can add numerous annotation tracks from the Tracks menu by choosing Annotation Tracks.

- - Ruler
 - Genes
 - RepeatMasker
 - Conservation
 - Genome Annotation

Each header can be expanded to one or more submenus that display tracks that can be added to the browser. The tracks include CpG island information, repeat information, G/C content information, and conservation information to name a few.

-	m	m10
		Ruler
		Genes
		RepeatMasker
		Conservation
	-	Genome Annotation
		mm10 Encode Blacklist
		CpG Context Add
		CpG Islands Add
		CpGs Add
		GC percent Add

1.6.3 Adding a custom track or data hub

Add

Users can also submit their own track as a custom track. For example, say we have a bigWig track located at https: //vizhub.wustl.edu/public/tmp/TW463_20-5-bonemarrow_MeDIP.bigWig . From the Tracks menu choose **Custom tracks** and a custom track interface will pop up. Fill in the track type, label, and URL before clicking the green *Submit* button:

Add Custom Track

Add Custom Data Hub

Add custom track

Track type

bigWig - numerical data

Track label

MeDIP

Track file URL

https://wangftp.wustl.edu/~dli/test/TW463_20-5-bonemarrow_MeDIP.bigWig

Submit

You can see the track is added:

Tools:	些	24	- +	Ø	C	<u>1</u>		
Ruler		.0M		52.1	M		52.2M	52.3M
MeDIP		10 0	and the state	a de se	An	A.	A A.	والمحدوم الأفرور أتشط والمسروح

Adding a custom data hub is similar to the steps above. For example, say you have a hub located at https://vizhub. wustl.edu/public/tmp/a.json . From the Tracks menu choose **Custom tracks**, switch to the *Add custom data hub* tab, paste the URL of your hub, and then click the green *Load From URL* button. from URL.

Add Custom Track	Add Custom Data Hub			
Add custo	m data hub			
Custom hub URL				
https://wangftp.wustl.ed	du/~dli/test/a.json			
Load from URL				
Or				
Choose datahub file				Browse
Row:	sample	ţ	Column: assay	
			massay	
	⊞sample		0/2	

The tracks within the custom hub can then be added from the generated facet table.

Note: Tracks from custom hubs are hidden by default as users may submit a hub contains hundreds of tracks. Users should add tracks that they want from the facet table.

You can also load a local data hub file in JSON format from your computer using the *file upload* interface, right below the *URL submit* hub interface.

Also see the Tracks and Datahub sections for how to prepare your tracks and datahub files.

1.7 Track Customization

Tracks can be customized in a multitude of manners.

1.7.1 Selecting Tracks

An indivdual track can be selected by simply right clicking on the tracking on the track. Multiple tracks can be selected by either holding the *shift* button and left clicking on each track or by holding shift and left clicking on a shared metadata term of consecutive tracks. In this manner, multiple tracks can be customized or moved at the same time. To deselect the tracks simply right click and press the button Deselect # tracks.

1.7.2 Track Color

Right clicking on annotation and numerical tracks will display Primary Color, Secondary Color, and Background Color which can all be customized using the color picker. For methylC tracks and categorical tracks the Color and Background of each class of elements (e.g. CG, CHG, and CHH) can be personalized. Additionally, for methylC``tracks the ``Read depth line color can be customized.

1.7.3 Track Height

For each track the height can be customized by right clicking on the track and typing in a number to the panel. At 20 pixels and below for numerical tracks the track will display as a heatmap.

1.7.4 Track Display Mode

For each numerical, annotation, or BAM track the display can be changed to DENSITY or FULL mode by right clicking on the track.

1.7.5 Track Y-axis Scale

For each numerical track the y-axis can be displayed in AUTO or FIXED mode by right clicking on the track. The AUTO mode will scale the axis based on numerical values in the immediate area of the view range. The FIXED mode allows the user to select a Y-Axis min or Y-axis max. For values above the set max the Primary color above max can be set for easy viewing. For values below the set minimum the Primary color below min can bet set.

1.7.6 Track Information

If details were specified for a track in the data hub file these can be viewed by right clicking on the sample and clicking on the arrow to the right. An easy access copy but is also available to copy the URL for the track.



1.8 Circlet view for chromatin interaction tracks


For any chromatin interaction track type (*HiC*, *longrange*, *bigInteract*), when you right click the track, you can see the Circlet view button:

¥ .			
	•	test	
		Display mode: HEATMAP •	
	•	Score scale: AUTO 🔻	
		Primary color	
	•	Secondary color	
•		Background color	
		Circlet view	1. 10 A
		× Remove	
		More information	
		•	•
			•

Click the button will bring you to the Circlet view interface. You can config the layout and/or the data source:

Choose a layout range:				
Current region C	Current chromosome	0	Whole genome	9
Choose data from: (please note only .hic track will fetch	additional data)			
Current region 9	Current chromosome	0	Whole genome	0
Change color: #a813c7				
Change score scale:				
Min: 1 Max	10			
Change flanking at each si	de:			
Length: 180000	bp			
Turn off scale 🗹				
1 Download				

And config the color, scale, flanking region length at each end of one interaction. You also can download the view as a SVG file used for publication.



CHAPTER 2

Tracks

Important: Since all tracks are hosted on the web with HTTP/HTTPS links provided for submission as tracks, the webservers which are hosting the track files need Cross-Origin Resource Sharing (CORS) enabled.

Quoted from MDN:

```
Cross-Origin Resource Sharing (CORS) is a mechanism that uses additional HTTP headers to tell a browser to let a web application running at one origin (domain) have permission to access selected resources from a server at a different origin. A web application makes a cross-origin HTTP request when it requests a resource that has a different origin (domain, protocol, and port) than its own origin.
```

2.1 Configure your webserver to enable CORS

Most likely the browser domain is different from the server the tracks are hosted on. The hosting server needs CORS enabled and for an Apache web server in Ubuntu this setup will work:

```
Header always set Access-Control-Allow-Origin "*"
Header always set Access-Control-Allow-Methods "POST, GET, OPTIONS, DELETE, PUT"
Header always set Access-Control-Max-Age "1000"
Header always set Access-Control-Allow-Headers "x-requested-with, Content-Type,

→origin, authorization, accept, client-security-token"
```

2.2 Prepare track files

The browser accesses track files from their URL. Only a portion of the data, that within the specific view region, are transferred to the browser for visualization. Thus, all the track files need be hosted in a web accessible location using HTTP or HTTPS. The following sections introduce the track types that the browser supports.

Binary track file formats like *bigWig* and *HiC* can be used directly with the browser.

bedGraph, methylC, categorical, longrange and *bed* track files need to be compressed by bgzip and indexed by tabix for use by the browser. The resulting index file with suffix .tbi needs to be located at the same URL with the .gz file.

Bed like format track files need be sorted before submission. For example, if we have a track file named track. bedgraph we can use the generic Linux sort command, the bedSort tool from UCSC, or the sort-bed command from BEDOPS. Here is an example command using each of the three methods:

```
# Using Linux sort
sort -k1,1 -k2,2n track.bedgraph > track.bedgraph.sorted
# Using bedSort
bedSort track.bedgraph track.bedgraph.sorted
# Using sort-bed
sort-bed track.bedgraph > track.bedgraph.sorted
```

Then the file must be compressed using bgzip and indexed using tabix:

```
bgzip track.bedgraph.sorted
tabix -p bed track.bedgraph.sorted.gz
```

Move files "track.bedgraph.sorted.gz" and "track.bedgraph.sorted.gz.tbi" to a web server. The two files must be in the same directory. Obtain the URL to "track.bedgraph.sorted.gz" for submission.

SAM files first need to be compressed to *BAM* files. *BAM* files need to be coordinate sorted and indexed for use by the browser. The resulting index file with suffix .bai needs be located at the same URL with the .bam file.

Here is an example command:

```
# Using samtools view to convert to bam
samtools view -Sb test.sam > test.bam
# Using samtools sort to coordinate sort the file
samtools sort test.bam > test.sorted.bam
# Using samtools index
samtools index test.sorted.bam
```

2.3 Annotation Tracks

Annotation tracks represent genomic features or intervals across the genome. Popular examples include SNP files, CpG Island files, and blacklisted regions.

2.3.1 bed

bed format files can be used to annotate elements across the genome or to represent reads from a sequencing experiment. For more about the bed format please check the UCSC bed page.

Example lines are below:

 chr9
 3035610
 3036180
 Blacklist_155
 +

 chr9
 3036200
 3036480
 Blacklist_156
 +

 chr9
 3036420
 3036660
 Blacklist_157
 +

Every line must consist of at least 3 fields separated by the Tab delimiter. The required fields from left to right are chromosome, start position (0-based), and end position (not included). A fourth (optional) column is

reserved for the name of the interval and the sixth column (optional) is reserved for the strand. All other columns are ignored, but can be present in the file.

Blacklist 3 Column	
Blacklist 4 Column	Blacklist 155
Blacklist 6 Column Positive	>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>
Blacklist 6 Column Negative	K<<<<<< Blacklist 155 <<<<<< K<<<<<<<<< K<<<<<<<<<<<<<<<<<<

Note: The display of a bed file differs by how many columns are provided in the file (see image above). The simplest, 3 column, format just displays blocks for each interval. The four column format displays the name of each element over each interval. If the sixth column is provided in the file then >>> or <<< will be displayed over each interval to represent strand information.

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

2.3.2 refbed

The refbed format files allows you to upload a custom gene annotation track. It is similar to the refGene bed-like file downloaded from UCSC but with slight modifications. Each file of this format contains (each column is separated by *Tab*):

chr, transcript_start, transcript_stop, translation_start, translation_stop, strand, gene_name, transcript_id, type, exon(including UTR bases) starts, exon(including UTR bases) stops, and additional gene info (*optional*)

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

Hint: The 9th column contains gene type, but is simplified from the Gencode/Ensembl annotations to coding, pseudo, nonCoding, problem, and other. These classes of gene type are colored differently when the track is displayed on the browser.

Hint: The 10th and 11th columns contain exon starts and ends respectively. Each start or end is seperated by a comma.

For example:

```
start1,start2,start3,start4 stop1,stop2,stop3,stop4
100,120,140,160 110,130,150,170
```

Hint: The 12th column contains extra information. This information can be manually annotated or we suggest using Ensembl Biomart to download paired Transcript stable IDs and Gene descriptions. The information in this column must be separated by *spaces* and not tabs.

All of the below lines will work for additional information in the 12th column:

Gene ID:ENSMUSG00000103482.1 Gene Type:TEC Transcript Type:TEC Additional →Info:predicted gene, 37999 [Source:MGI Symbol;Acc:MGI:5611227] Gene ID:ENSMUSG00000103482.1 Gene Type:TEC Transcript Type:TEC ENSMUSG00000103482.1 TEC Additional Info:predicted gene, 37999 [Source:MGI Symbol;Acc:MGI:5611227] My Favorite Gene

Here are a few example lines in refbed format from gencode.vM17.annotation.gtf (mouse mm10 format):

```
24910461
                          24911659
                                         24910461
                                                        24911659
chr1
                                       pseudo 24911220,24910461
→RP23-109H7.1
               ENSMUST00000187022.1
                                                                     24911659,
                         ID:ENSMUSG00000100808.1 Gene Type:processed_pseudogene
→24910681
               Gene
→Transcript Type:processed_pseudogene Additional Info:predicted gene 28594
→ [Source:MGI Symbol; Acc:MGI:5579300]
         25203443
                                         25203443
                                                        25205696
chr1
                          25205696
                                                                       _
→Adgrb3 ENSMUST00000190202.1
                               coding 25203443
                                                      25205696
                                                                     Gene
                     ID:ENSMUSG00000033569.17 Gene Type:protein_coding Transcript_
→Type:retained_intron Additional Info:adhesion G protein-coupled receptor
                                                                          ВЗ.
→ [Source:MGI Symbol; Acc:MGI:2441837]
                                         25276404
chr1
         25276404
                         25277954
                                                        25277954
                                                                              problem 25276404
                                                             25277954
               ENSMUST00000193138.1
→RP23-21P2.4
                            ID:ENSMUSG00000104257.1 Gene Type:TEC Transcript
⊶Gene
→Type:TEC Additional Info:predicted gene, 20172 [Source:MGI Symbol;Acc:MGI:5012357]
                                        26566833 26566938
chr1 26566833 26566938
                                                                       +
→Gm24064 ENSMUST00000157486.1
                                            26566833
                                                             26566938
                              nonCoding
→Gene
                              ID:ENSMUSG0000088111.1 Gene Type:snoRNA Transcript,
→Type:snoRNA Additional Info:predicted gene, 24064 [Source:MGI
→ Symbol;Acc:MGI:5453841]
```

Note: The last optional column is dislayed as a gene description when a gene is clicked on the browser. Our modified format can be easily obtained from available refGene.bed file downloads from UCSC. Gencode GTF and Ensembl GTF files can be manipulated to this format using the Converting_Gencode_or_Ensembl_GTF_to_refBed.bash script in scripts. The script by default puts Gene ID:, Gene Type:, and Transcript Type: in the additional information column. Run with an annotation file, with columns Transcript_ID and Description (seperated by a tab), the script will also add "Additional Info" to the 12th column. The script depends on bedtools, bgzip, and tabix. Lastly, within the script an awk array is used to reclassify gene type and can easily be modified for additional gene types.

The script is run as follows:

```
bash Converting_Gencode_or_Ensembl_GTF_to_refBed.bash Ensembl my.gtf my_optional_

→annotation.txt
bash Converting_Gencode_or_Ensembl_GTF_to_refBed.bash Gencode gencode.vM17.annotation.

→gtf
bash Converting_Gencode_or_Ensembl_GTF_to_refBed.bash Gencode gencode.vM17.annotation.

→gtf biomart_2col.txt
```

Warning: Spaces are used as delimiters in the GTF files so change gene names with spaces before processing.

For Example:

sed -i 's/ (1 of many)/_(1_of_many)/g' Danio_rerio.GRCz10.91.chr.gtf

2.4 Numerical Tracks

Currently there are two types of numerical tracks:

- bigWig
- bedGraph

2.4.1 bigWig

bigWig is a popular format to represent numerical values over genomic coordinates. Please check the UCSC bigWig page to learn more about this format.

2.4.2 bedGraph

bedGraph format also defines values in different genomic locations. For more about the bedGraph format please check the UCSC bedGraph page.

Example lines are below:

```
chr12 6537598 6537599 28.80914
chr12 6537599 6537600 28.96908
chr12 6537599 6537612 -2
chr12 6537600 6537601 29.30229
```

Every line consists of 4 fields separated by the Tab delimiter. The fields from left to right are chromosome, start position (0-based), end position (not included), and value.

Note: You can use negative values for reverse strand. Both positive and negative values can exist over the same coordinates (they can overlap). In bigWig format negative values can also be specified, but they cannot overlap with positive values.

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

2.5 Read Alignment BAM Tracks

2.5.1 BAM

The BAM format is a compressed SAM format used to store sequence alignment data. Please check the Samtools Documentation page to learn more about this format and how to manipulate these files.

2.6 Methylation tracks

Methylation experiments like MeDIP-seq or MRE-seq can use *bigWig* or *bedGraph* format for data display. For WGBS if users want to show read depth, methylation context, and methylation level then the data is best suited for the *methylC* format, described below.

2.6.1 methylC

chr1	10542	10543	CG	0.923	-	26
chr1	10556	10557	CHH	0.040	_	25
chr1	10562	10563	CG	0.941	+	17
chr1	10563	10564	CG	0.958	_	24
chr1	10564	10565	CHG	0.056	+	18
chr1	10566	10567	CHG	0.045	_	22
chr1	10570	10571	CG	0.870	+	23
chrl	10571	10572	CG	0.913	-	23

Methylation data are formatted in methylC format, which is a 7 column bed format file:

Each line contains 7 fields separated by Tab. The fields are chromosome, start position (0-based), end position (not included), methylation context (CG, CHG, CHG etc.), methylation value, strand, and read depth.

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

2.7 Categorical tracks

Categorical tracks represent genomic bins for different categories. The most popular example is the representation of chromHMM data which indicates which region is likely an enhancer, likely a promoter, etc. Other uses for the track include the display of different types of methylation (DMRs, DMVs, LMRs, UMRs, etc.) or even peaks colored by tissue type.

2.7.1 categorical

The categorical track uses the first three columns of the standard *bed* format (chromosome, start position (0-based), and end position (not included)) with the addition of a 4th column indicating the category type which can be a string or number:

chr1	start1	end1	category1
chr2	start2	end2	category2
chr3	start3	end3	category3
chr4	start4	end4	category4

Important: when you use numbers like 1, 2 and 3 as category names, in the datahub definition, please use it a string for the category attribute in options, see the example below:

```
{
    "type": "categorical",
    "name": "ChromHMM",
    "url": "https://egg.wustl.edu/d/hg19/E017_15_coreMarks_dense.gz",
    "options": {
        "category": {
            "1": {"name": "Active TSS", "color": "#ff0000"},
            "2": {"name": "Flanking Active TSS", "color": "#ff4500"},
            "3": {"name": "Transcr at gene 5' and 3'", "color": "#32cd32"}
    }
}
```

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

2.8 Long range chromatin interaction

Long range chromatin interaction data are used to show relationships between genomic regions. HiC is used to show the results from a HiC experiment.

2.8.1 HiC

To learn more about the HiC format please check https://github.com/aidenlab/juicer/wiki/Data.

2.8.2 longrange

The longrange track is a *bed* format-like file type. Each row contains columns from left to right: chromosome, start position (0-based), and end position (not included), interaction target in this format chr2:333-444, 55. As an example, interval "chr1:111-222" interacts with interval "chr2:333-444" on a score of 55, we will use following two lines to represent this interaction:

chr1 111 222 chr2:333-444,55 chr2 333 444 chr1:111-222,55

Important: Be sure to make **TWO** records for a pair of interacting loci, one record for each locus.

This format needs to be compressed by bgzip and indexed by tabix for submission as a track. See Prepare track files.

2.8.3 bigInteract

The bigInteract format from UCSC can also be used at the browser, for more details about this format, please check the UCSC bigInteract format page.

2.8.4 cool

Thanks to the higlass team who provides the data API, the browser is able to display cool tracks by using the data uuid from the higlass server, for example, you can use the uuid Hyc3TZevQVm3FcTAZShLQg to represent the track for *Aiden et al. (2009) GM06900 HINDIII 1kb*, for a full list of available cool tracks please check http://higlass.io/api/v1/tilesets/?dt=matrix

CHAPTER 3

Datahub

A datahub is a JSON file descibing a set of tracks in the browser. A datahub file is an array of tracks, which are also defined in JSON syntax:

```
[
    "type": "track_type1",
    "name": "track_name1",
    "url": "track_url1",
    "options": {
        "color": "red"
        }
    },
    {
    "type": "track_type2",
    "name": "track_name2",
    "url": "track_url2",
    "options": {
        "color": "blue"
        }
    }
]
```

3.1 Example data hub

Pasted below is an example data hub for mouse mm10:

```
[
    {
        "type": "bigwig",
        "url": "https://vizhub.wustl.edu/public/tmp/TW463_20-5-bonemarrow_MeDIP.bigWig",
        "name": "MeDIP",
        "options": {
        "options": {
        }
    }
}
```

(continues on next page)

(continued from previous page)

```
"color": "red",
        "backgroundColor": "#FFE7AB"
        },
    "metadata": {
        "sample": "bone",
        "assay": "MeDIP"
        }
    },
    {
    "type": "bigwig",
    "url": "https://vizhub.wustl.edu/public/tmp/TW551_20-5-bonemarrow_MRE.CpG.bigWig",
    "name": "MRE",
    "options": {
        "color": "blue",
        "backgroundColor": "#C0E3CC"
        },
    "metadata": {
        "sample": "bone",
        "assay": "MRE"
        }
    }
1
```

3.2 Example bigWig track

```
{
   "type": "bigwig",
   "name": "example bigwig",
   "url": "https://vizhub.wustl.edu/hubSample/hg19/GSM429321.bigWig",
   "options": {
        "color": "blue"
    }
}
```

3.3 Example methylC track

```
{
    "type": "methylc",
    "name": "H1",
    "url": "https://vizhub.wustl.edu/public/hg19/methylc2/h1.liftedtohg19.gz",
    "options": {
        "label": "Methylation",
        "colorsForContext": {
            "CG": { "color": "#648bd8", "background": "#d9d9d9" },
            "CHG": { "color": "#ff944d", "background": "#ffe0cc" },
            "CHH": { "color": "#ff00ff", "background": "#ffe5ff" }
        },
        "depthColor": "#01E9FE"
    },
}
```

3.4 Example categorical track

```
"type": "categorical",
 "name": "ChromHMM",
 "url": "https://egg.wustl.edu/d/hg19/E017_15_coreMarks_dense.gz",
 "options": {
     "category": {
         "1": {"name": "Active TSS", "color": "#ff0000"},
         "2": {"name": "Flanking Active TSS", "color": "#ff4500"},
         "3": {"name": "Transcr at gene 5' and 3'", "color": "#32cd32"},
         "4": {"name": "Strong transcription", "color": "#008000"},
         "5": {"name": "Weak transcription", "color": "#006400"},
         "6": {"name": "Genic enhancers", "color": "#c2e105"},
         "7": {"name": "Enhancers", "color": "#ffff00"},
         "8": {"name": "ZNF genes & repeats", "color": "#66cdaa"},
         "9": {"name": "Heterochromatin", "color": "#8
                                                         a91d0"},
         "10": {"name": "Bivalent/Poised TSS", "color": "#cd5c5c"},
         "11": {"name": "Flanking Bivalent TSS/Enh", "color": "#e9967a"},
         "12": {"name": "Bivalent Enhancer", "color": "#bdb76b"},
         "13": {"name": "Repressed PolyComb", "color": "#808080"},
         "14": {"name": "Weak Repressed PolyComb", "color": "#c0c0c0"},
         "15": {"name": "Quiescent/Low", "color": "#ffffff"}
     }
 }
}
```

Supported options: backgroundColor, color, color2, yScale, yMax, and yMin.

3.5 Example longrange track

```
"type": "longrange",
"name": "ES-E14 ChIA-PET",
"url": "https://egg.wustl.edu/d/mm9/GSE28247_st3c.gz"
```

3.6 Example bigInteract track

```
{
    "type": "biginteract",
    "name": "test bigInteract",
    "url": "https://epgg-test.wustl.edu/dli/long-range-test/interactExample3.inter.bb"
}
```

3.7 Example repeatmasker track

```
"type": "repeatmasker",
"name": "RepeatMasker",
"url": "https://vizhub.wustl.edu/public/mm10/rmsk16.bb"
```

3.8 Example geneAnnotation track

```
"type": "geneAnnotation",
"name": "refGene",
"genome": "mm10"
```

{

}

{

}

Note: Please specify the genome attibute for gene annotation tracks.

3.9 Example bigbed track

```
"type": "bigbed",
"name": "test bigbed",
"url": "https://vizhub.wustl.edu/hubSample/hg19/bigBed1"
```

3.10 Example bed track

```
{
    "type": "bed",
    "name": "mm10 bed",
    "url": "https://epgg-test.wustl.edu/d/mm10/mm10_cpgIslands.bed.gz"
```

3.11 Example refbed track

```
"type": "refbed",
"name": "mm10 gencode basic",
"url": "https://vizhub.wustl.edu/public/tmp/gencodeM18_load_basic_Gene.bed.gz",
"options": {
        "categoryColors": {
            "coding": "rgb(101,1,168)",
            "nonCoding": "rgb(1,193,75)",
            "pseudo": "rgb(230,0,172)",
            "problem": "rgb(224,2,2)",
```

(continues on next page)

{

(continued from previous page)

```
"other":"rgb(128,128,128)"
}
}
```

Note: categoryColors designates colors for the gene type as indicated in the 9th column. The default scheme is shown above for the five classes created by the Converting_Gencode_or_Ensembl_GTF_to_refBed. bash script, but any number of categories can be defined.

3.12 Example HiC track

{

}

```
"type": "hic",
"name": "test hic",
"url": "https://epgg-test.wustl.edu/dli/long-range-test/test.hic",
"options": {
    "displayMode": "arc"
}
```

3.13 Example cool track

```
{
  "type": "cool",
  "name": "Aiden et al. (2009) GM06900 HINDIII 1kb",
  "url": "Hyc3TZevQVm3FcTAZShLQg",
  "options": {
     "displayMode": "arc"
  }
}
```

Note: please note we are using the uuid Hyc3TZevQVm3FcTAZShLQg here from higlass API server instead of a file URL to represent a cool track.

3.14 Example genomealign track

```
{
   "name": "hg19 to mm10 alignment",
   "type": "genomealign",
   "metadata": {
        "genome": "mm10"
    }
}
```

3.15 Example Ruler track

```
"type": "ruler",
"name": "Ruler"
```

3.16 Track properties

3.16.1 type

{

}

Requiied. type specifies the track type, currently supported track types:

- bigWig
- bedGraph
- methylC
- categorical
- hic
- bed
- bigbed
- refbed
- repeatmasker
- geneAnnotation
- genomealign
- longrange
- bigInteract
- ruler

Note: type is case insensitive.

3.16.2 name

Required. name specifies the track name used internally by the browser. It is also displayed as the track legend if no *label* specificied. Value can be any string.

3.16.3 label

Optional. label specifies the track legend displayed in the browser. It overrides the *name* artibute. Value can be any string.

3.16.4 url

Required. url contains the URL to the track file and needs to be HTTP or HTTPS location string.

```
Important: A url is required for all the tracks in binary format. Gene annotaion tracks, like refGene, do not need a url as they are stored in the Mongo database. Additional annotation tracks, such as the ruler track, also do not need a url.
```

Caution: Each user-provided url must link to a publically available website, without password protection, so that the browser can read in the file.

3.16.5 metadata

Optional. An object specifying the metadata of the track.

In this basic example the value of each metadata term is a string.

```
"metadata": {
    "sample": "bone",
    "assay": "MRE"
}
```

This example public Roadmap data hub has more complex metadata definitions and makes use of a **list of strings** to build a *hierarchical structure*.

```
{
    "url": "https://egg.wustl.edu/d/hg19/GSM997242_1.bigWig",
    "metadata": {
        "Sample": [
            "Adult Cells/Tissues",
            "Blood",
            "Other blood cells",
            "CD4+_CD25-_Th_Primary_Cells"
        ],
        "Donor": [
            "Donor Identifier",
            "Donor_332"
        ],
        "Assay": [
            "Epigenetic Mark",
            "Histone Mark",
            "нз",
            "H3K9",
            "H3K9me3"
        ],
        "Institution": [
            "Broad Institute"
        1
   },
    "type": "bigwig",
    "options": {
        "color": "rgb(159,0,72)"
   },
```

(continues on next page)

(continued from previous page)

```
"name": "H3K9me3 of CD4+_CD25-_Th_Primary_Cells"
```

The list of metadata is ordered from more generic to more specific and helps build the facet table hierarchy making the **search** and **filter** functions in track table easier.

3.16.6 details

Optional. If you want to add more information for each track then the details attribute is helpful. After right clicking on the track you can click **More Information** and see the details, url, and metadata for each track in the dropdown menu.

3.16.7 options

Optional. All track render options are placed in an object called options. This object can have the following properties:

color

color is used to define the color for each track. A color name, RGB values, or hex color code can be used. For more about color name or RGB please see https://www.w3schools.com/css/css_colors.asp.

color2

color2 is used to define the color for negative values from the track data. The default is the same as *color*.

backgroundColor

backgroundColor defines the background color of the track.

height

height controls the height of the track which is specified as a number and displayed in pixels.

yScale

yScale allows you to configure the track's y-scale. Options include *auto* or *fixed. auto* sets the y-scale from 0 to the max value of values in the view region for a given track. *fixed* means you can specify the *minimal* and *maximal* value.

yMax

yMax is used to define the maximum value of a track's y-axis. Value is number.

yMin

yMin is used to define the *minimum* value of a track's y-axis. Value is number.

Important: If you need the track to be in *fixed* scale, you need to specify yScale to *fixed* besides of set yMax and yMin.

colorAboveMax

colorAboveMax defines the color displayed when a *fixed yScale* is used and a value exceeds the *yMax* defined.

color2BelowMin

color2BelowMin defines the color displayed when a *fixed yScale* is used and a value is below the *yMin* defined.

displayMode

displayMode specifies display mode for each tracks. Different tracks have different display modes as listed below.

type	display mode
bigWig	auto, bar, heatmap
bedGraph	auto, bar, heatmap
geneAnnotation	full, density
HiC	arc, heatmap
genomealign	rough, fine

Installation

4.1 Setup

• Install NodeJS from https://nodejs.org/en/

Note: Feel free to use any package manager tool on your system for installation (brew, etc.).

4.2 Start the browser

- 1. Enter the frontend directory
- 2. Type npm install (just for the first time) This step will install dependent packages.
- 3. Type npm start

That's it! You are done with your mirror site. The browser is now accessible from http://localhost:3000/browser.

4.3 Setup your own backend API (optional)

By default, your local browser mirror site uses our API service at https://api.epigenomegateway.org/ documentation, while if you find the species or assembly you are interested is not listed by our API, you can either contact us to add it or build your own API. To build your own API, please follow the steps below:

- 1. Install MongoDB from https://www.mongodb.com/
- 2. start mongdb server
- 3. Enter the backend directory
- 4. Type npm install

Then prepare your gene annotation files like the ones for hg19, mm10 etc:

- 1. Make sure MongoDB is running
- 2. Enter the backend directory
- 3. Run npm run setup This step will load the gene annotation data to the MongoDB database
- 4. Type npm start

Now your own backend API is running, change AWS_API variable to empty string in GeneSource.js file. After this you are using your own API for gene annotation tracks and gene search.

Embedding

To embed the browser in any HTML file, create a HTML page with following contents: (the example shows how to embed a mouse browser with 2 bigWig tracks from ENCODE data portal)

```
<html lang="en">
<head>
<meta charset="utf-8">
<meta name="viewport" content="width=device-width, initial-scale=1, shrink-to-fit=no">
<meta name="theme-color" content="#000000">
<title>The New WashU Epigenome Browser</title>
<link rel="stylesheet" href="https://maxcdn.bootstrapcdn.com/bootstrap/4.0.0/css/</pre>
→bootstrap.min.css" integrity="sha384-
→Gn5384xqQ1aoWXA+058RXPxPq6fy4IWvTNh0E263XmFcJlSAwiGqFAW/dAiS6JXm"
    crossorigin="anonymous">
<script src="https://iqv.org/web/release/2.0.1/dist/iqv.min.js"></script>
<script src="https://iqv.org/web/jb/release/1.0.0/dist/juicebox.min.js"></script>
<script src="https://aframe.io/releases/0.8.0/aframe.min.js"></script>
<script src="https://code.jquery.com/jquery-3.2.1.slim.min.js" integrity="sha384-
→KJ3o2DKtIkvYIK3UENzmM7KCkRr/rE9/Qpg6aAZGJwFDMVNA/GpGFF93hXpG5KkN"
   crossorigin="anonymous"></script>
script src="https://cdnjs.cloudflare.com/ajax/libs/popper.js/1.12.9/umd/popper.min.js
→ " integrity="sha384-ApNbgh9B+Y1QKtv3Rn7W3mgPxhU9K/ScQsAP7hUibX39j7fakFPskvXusvfa0b4Q
\hookrightarrow "
    crossorigin="anonymous"></script>
<script src="https://maxcdn.bootstrapcdn.com/bootstrap/4.0.0/js/bootstrap.min.js"...</pre>
→integrity="sha384-JZR6Spejh4U02d8jOt6vLEHfe/JQGiRRSQQxSfFWpi1MquVdAyjUar5+76PVCmY1"
   crossorigin="anonymous"></script>
<script src="https://unpkg.com/epgg@latest/umd/epgg.js"></script>
<link rel="stylesheet" href="https://unpkg.com/epgg@latest/umd/epgg.css">
</head>
<body>
<noscript>
   You need to enable JavaScript to run this app.
</noscript>
<h1>Embedding test</h1>
```

(continues on next page)

(continued from previous page)

```
<div id="embed" style="width:1000px"></div>
<h2>some other headings</h2>
<script>
   const container = document.getElementById('embed');
    const contents = {
       "genomeName": "mm10",
        "displayRegion": "chr5:51997494-52853744",
        "trackLegendWidth": 120,
        "isShowingNavigator": true,
        "tracks": [
        {
            "type": "geneannotation",
            "name": "refGene",
            "genome": "mm10"
        },
        {
            "type": "geneannotation",
            "name": "gencodeM19Basic",
            "genome": "mm10"
        },
        {
            "type": "ruler",
            "name": "Ruler"
        },
        {
            "type": "bigWig",
            "name": "ChipSeg of Heart",
            "url": "https://www.encodeproject.org/files/ENCFF641FBI/@@download/
→ENCFF641FBI.bigWig",
            "options": { "color": "red" },
            "metadata": { "Sample": "Heart" }
        },
        {
            "type": "bigWig",
            "name": "ChipSeq of Liver",
            "url": "https://www.encodeproject.org/files/ENCFF555LBI/@@download/
→ENCFF555LBI.bigWig",
           "options": { "color": "blue" },
            "metadata": { "Sample": "Liver" }
        }
       ],
        "metadataTerms": ["Sample"],
        "regionSets": [],
        "regionSetViewIndex": -1,
    };
    renderBrowserInElement(contents, container);
</script>
</body>
</html>
```

The key API is the function renderBrowserInElement, it accepts the contents array as first argument, and container as second argument which is a DOM element.

Frontend code architeture

Note: This section explains how frontend code is organized, intend to be used for development purpose. Regular browser users don't need to care about this section.

6.1 Quick tour

The client code is in the frontend folder. Here is a quick tour of frontend/src:

- components: All React components.
 - genomeNavigator: The navigation bar at the top that allows users to navigate
 - track: Track-related components
 - trackManagers: UI that manages adding tracks
- dataSources: API calls, AJAX calls, database connections, etc. that get data to display.
- model: Data models.
- stories: Stories for Storybook on which unit tests depend.
- vendor: 3rd-party libraries that are not in NPM.

6.2 Suggested order of reading

If you plan to understand the app as a whole here is a suggested order to read the code in:

- 1. Feature: A feature or annotation in the genome.
- 2. NavigationContext: A list of Features that represent everywhere a user can navigate. If the Features are actually entire chromosomes then the user can effectively navigate the whole genome.

- 3. DisplayedRegionModel: An interval in a NavigationContext.
- 4. App: The root component of the app.
- 5. From App, descend into interested components.

6.3 Making a new track type

6.3.1 Make a new TrackConfig

Make a new class that extends TrackConfigor one of its subclasses. This class packages many essential track characteristics:

- getComponent() Gets the component that renders the main visualizer and legend of the track.
- getMenuComponents () Specifies context menu items in an array of components. You can choose existing ones in the contextMenu directory or make new ones.
- getOptions() The visualizer probably renders with default options like a color. This method returns a plain object containing those options.

You do not have to implement these methods immediately as the base TrackConfig class provides minimal defaults. Just work on making the browser render *some* temporary placeholder at first.

6.3.2 Specify when to use the TrackConfig

- 1. Import your new TrackConfig into trackConfig/getTrackConfig.js.
- 2. Add an appropriate entry to TYPE_NAME_TO_SUBTYPE, which maps track type name to track renderer.

6.3.3 Write a new track visualizer component (implement getComponent())

- 1. Make a new component expecting to receive a bunch of props from TrackContainer. Track.js documents the props to expect.
- 2. If you need data assume it will come through the data prop. We will add data fetch in the next step.
- 3. Your new component may render anything though it is **highly** recommended you render a <Track> component, if not one of the more specialized components like <AnnotationTrack> or <NumericalTrack>. Pass *all* track container props to these sub-components.
- 4. In addition to track container props you need to provide certain props to these sub-components, all of which the respective files document.
 - For example, <Track> requires a legend and visualizer element. Use the track container props, which includes view region and width, to render a visualizer and pass it to <Track>.

6.3.4 Add data fetch

Available data sources are in the dataSources folder. If none of them fulfill your needs, write a new class that fulfills the interface of DataSource.js. More can be found in that file.

How do we give your visualizer data? Higher-order components! track/commonComponents contains track-specific HOCs; their names start with config- or with-.

configStaticDataSource requests a callback that returns a DataSource and then returns a *function* that wraps React components. After you use this function, a component will automatically receive three props data, isLoading, and error. These update with the browser's current view region. In particular, the HOC guarantees synchronization of the data prop with the current view region if isLoading is false.

6.3.5 2. Specify context menu components (implement getMenuComponents())

Specify context menu items with an array of components. You can choose existing ones in the contextMenu directory or make new ones.

• Make sure the method returns Component classes, not component instances.

6.3.6 3. Specify default options

Default option objects look like the options prop of TrackModel objects. Context menu items will read these options if the track model does not specify them. Make sure these options are consistent with the way you are rendering your track component! The configOptionMerging HOC should help with that.

Once you have a default options object, call setDefaultOptions () in the constructor of TrackConfig to use them.

6.4 Performance tips

Querying the width or height of any element, for example through clientWidth or getBoundingClientRect() is slow. Such queries take on the order of 2 to 20 ms. While it is fine to do it once or twice, avoid doing it in a loop. Suppose you aim to plot 500 data points on a SVG and for each point you query the SVG's width. That is already a second or more of computation – very noticable to the user!

6.5 React (and other) gotchas

- On Macs, control + click is the same as a right click which fires a contextmenu event. Note that click events do not fire on contextmenu events. The mousedown and mouseup events will still fire though.
- When using native DOM events they take priority over React events. This is because React waits for events to bubble to the root component before handling them. This can cause undesirable effects: for example, calling stopPropagation() on a React event will not actually stop native events. This StackOverflow post may also help if you have propagation problems: https://stackoverflow.com/questions/24415631/ reactjs-syntheticevent-stoppropagation-only-works-with-react-events
- React *always* unmounts components if their parents change type. The Reparentable component works around this by using app-unique IDs, but it can cause side effects with React's native events. Use with care.
- Webpack does not support circular dependencies, and while compilation may be successful, an import may resolve as undefined at runtime.

6.6 Lessons trying to refactor into WebWorkers

1. Data fetch and track display options are intimately related. For example, what if someone wants HiC data and selects the 5KB resolution option?

- 2. Thus, for each track type, we have one object that gets the track component, default rendering options, and data fetch/processing.
- 3. Webpack hangs forever if it encounters a cyclic dependency involving a webworker.
- 4. The code as in (2) causes a cyclic depdendency. This cycle is [config object] -> [data source] -> [worker] -> [track config deserializer] -> [config object]
- 5. We cannot have our cake and eat it too.

Unfortunately, this means we cannot pipeline all expensive computation in worker context, while also ensuring track component and data source live in the same place.

Contact us

7.1 Source code and issue tracker

You can find source code at our github repo: https://github.com/lidaof/eg-react

If you find any issue or have any question, please submit an issue here: https://github.com/lidaof/eg-react/issues. Thank you.

7.2 Get help on social media

- 1. Google plus forum: https://goo.gl/ruXs2Y
- 2. Facebook page: https://www.facebook.com/WashUEpiGenomeBrowser/
- 3. Twitter: https://twitter.com/wuepgg

CHAPTER $\mathbf{8}$

Indices and tables

- genindex
- search